# VI Package Manager

VI Package Manager (VIPM) is a community and tool that helps you find, download, share, & manage LabVIEW Add-ons.

# Resources

- https://vipm.io
- https://www.vipm.io/desktop/versions/
- https://forums.vipm.io

# FAQ

#### Package Cache

The downloaded VI packages can be found at

<ProgramData>\JKI\VIPM\cache

#### **Custom Actions**

You can define VIs to be called during install/uninstall, so-called custom actions. I've tried to copy zip files to C:\ProgramData\Hampel Software Engineering\PPLs instead of installing them into vi.lib and it turned out, that the Post-Install Custom Action is executed before the files were extracted!

### Linux

### Installation

Download the VIPM release for Linux: https://www.vipm.io/desktop/

Create two folders

- \$ sudo mkdir -p /usr/local/JKI/VIPM
- \$ sudo mkdir /etc/JKI

Change the owner of the folders to your current user, to gain write access

\$ sudo chown \$USER -R /usr/local/JKI/VIPM

\$ sudo chown \$USER /etc/JKI

Unzip to the new folder

\$ unzip vipm-22.1.2354-linux.zip -d /usr/local/JKI/VIPM/

### Usage

VIPM needs write permission to the LabVIW installation. Because LabVIEW is usually installed under /usr/local/natinst/ the normal user has no write access and VIPM needs to be executed as user "root". To do this on modern linux systems, it is remommended to use the command sudo or pkexec. Examples

\$ sudo /usr/local/JKI/VIPM/vipm

If this command returns an error "Unable to open x display", it is necessary to provide the environment variables to access the x display.

\$ pkexec env DISPLAY=\$DISPLAY XAUTHORITY=\$XAUTHORITY
/usr/local/JKI/VIPM/vipm

To start vipm with sudo vipm, you can create a symbolic link with this command

\$ sudo ln -s /usr/local/JKI/VIPM/vipm /usr/bin/vipm

*Alternative*: Change the owner of the LabVIEW folder to your current user, to gain write access and don't need sudo.

\$ sudo chown \$USER -R /usr/local/natinst/

From: https://dokuwiki.hampel-soft.com/ - HAMPEL SOFTWARE ENGINEERING

Permanent link: https://dokuwiki.hampel-soft.com/kb/labview/tools/vipm?rev=1692429110



Last update: 2023/08/19 07:11